

14T MRSI – protocol (created 20072021_061221_09122021)

date:	SCAN	GROUP	BDL
-------	------	-------	-----

study name

weight: 1708

study	TAIL_PRONE
	1H-Qsurf_Rat (phase 90)
	MRI_TxSuc
Time of start	12:26

FROM: MRSI_Coronal_09122021**1 Localizer 1** Adj: wobble (setup, stop)

E4 → Acquire localizer for position with **basic frequency/MR scan in the instruction cart**
fix the rat

2 Localizer 2 Adj: ref power (Gauss 6 kHz, 2mm slice top of the brain)

power= 0,0225 W

E6 → Acquire with **automatic instructions** → to have a shim before the next power adj

*Target power for MRSI:

New coil conf (phase 0) - 0.05-0.06 W
Coil config trad (phase 90) = 0.02-0.03 W

Loc MS **3_Localizer multislice_10_short** – for voxel and ellipsoid

E8 Adj: Ref power (Gauss 6 kHz, 2mm slice top of the brain)
Seq: 1 avg, 10 slices, Image size: 256x256, **FOV: 24x24** → position the slices
→ Acquire with **automatic instructions**

power: 0,0226 W

T2 turbo rare **4_T2 turbo rare for voxel position AXIAL (T2 turbo rare_6_54K)**

Seq: 2 avg, DS=4, Rare factor: 6, encoding start -1 (in resolution/encoding), 256x256, 20 slices, 0 gap, FOV=24x24
Position the slices, current shim

No = **E10** → Automatic acquisition**5_T2 turbo rare for voxel position CORONAL (T2 turbo rare_6_54K)**

Adj: B0 map (settings are already saved)

B0 map No: E900001

Seq: 2 avg, DS=4, Rare factor: 6, encoding start -1 (in resolution/encoding), 256x256, 20 slices, 0 gap

Position the slices +FOV = 24x24 as MRSI

Mapshim – elliptical 13.8x8x15, no iterative corr

No = **E12** → Automatic acquisitionCreate a mapshim report **SD values: 105,3 Hz → 33,3 Hz****STEAM** **6_STEAM_highres_JM_1109202_water_MRSI_01022021 (1st one)**

Adj: power calibration at the level of hippo (voxel)

position voxel - 10 x 2 x 10 mm

Seq: 16 avg, 1 rep, offset 0, no WS, 2 DS, TE=3, TM=10, TR=4

OVS (15/12mm, gap=0.3mm), spoilers 15-25-35%, ref scan 16, Seq spoilers : 25-30-20%

Current shim

No = **E15** → Automatic acquisition →topspin/tick in Reconstr. **lw = 24,589 Hz**

Save shim (no matter the lw always save the first shim)**7_STEAM_highres_JM_1109202_water_MRSI_01022021 – for shimming**

Adj: B0 map (settings are already saved)

B0 map No: E900002

Map shim on the voxel: cuboid, no margin, iterative corrections

No = **E17** → Automatic acquisition **lw = 24,589 → 28,324 Hz**

No = → If necessary repeat acq: B0 map and map shim

Final lw = 24,589 Hz

*Option: loc. Freq and loc shim in the ins. Cart - loc. Shim doesn't work for a big vox

Save shim**Optional: 8_STEAM_highres_JM_1109202_water_MRSI_01022021**

1 repetition with 32 averages, ref scan 16, OVS on

WS ON : WS pulses 84/150

→ If you change the WS note it here:

TE= 3ms

TM=10 ms

TR= 4s

No = → Automatic acquisition + Current shim

FID-MRSI**9_CSI_fid_FOVsat_WS_coronal_07122020**

Copy Slice orientation from Coronal RARE (scan 5)

Adjust the slice offset - Middle coordinate of the VOI Position

Adjust the position of saturation bands –**always perpendicular**

Use both Axial and Coronal RARE acq as support

FOV = 24x24, flip angle = 52, Dummy scans = 8

Linear encoding, start at -1, -1

weighted acquisition / standard acquisition – Always standard with 1 avg.

Preparation/WS - VAPORLaunch in the **setup mode** to find optimal WS and BW of the pulses (test also with 8x8 matrix)WS flip angles = **84/150**BW = **660 Hz**

*For the flip angles 84/150 and BW=660 works best so start with those

Number of averages: **1****Number of repetitions:** **1**No= **E20** → Automatic acquisition + current shim - **metabolites**

Duplicate sequence - preparation/WS - none

No= **E21** → Automatic acquisition + current shim - **water**

Matrix size= 31 x 31

Coronal Slice position: **8-9****E22: Linewidth check with 6-STEAM: lw = 23, 242 Hz**

If you repeat with other parameters

1st repeat**FID - MRSI**

weighted acquisition / standard acquisition - circle the one you select

Preparation/WS - VAPOR

Averages: **2**Repetitions: **1**No= **E23** → Automatic acquisition + current shim - **metabolites**No= **E24** → Automatic acquisition + current shim - **water****E25: Linewidth check: lw = 25, 146 Hz****2nd repeat****FID - MRSI**weighted acquisition / standard acquisition - circle the one you select **(After 2 hours)**

Preparation/WS - VAPOR

Averages = **1**Repetitions= **1**No= **26** → Automatic acquisition + current shim - **metabolites**No= **27** → Automatic acquisition + current shim - **water****E29: Linewidth check: lw = 25, 659 Hz****CENTER THE MRSI SLICE ON THE CORONAL IMAGE****E30: 2 avg - Met after 2 hours****E31: 2 avg - Wat after 2 hours****E32: Linewidth check: lw = 26, 589 Hz****E33: 1 avg - Met after 4 hours****E34: 1 avg - Wat after 4 hours**

