

28 02 2025 - Rat MRSI - High Res And 3D BA
SHAM - 793

9T MRSI - protocol created in PV360v3.5 based on protocols created on 14T: MRSI_2

date:	SCAN	GROUP	BDL
-------	------	-------	-----

study name **MAR-003067**

weight: **262,78**

INSIDE THE MAGNET TAPE THE HOLDER AND PUT THE BLUE LOCKER UPP

study	Application: MRSI_2_cryo/HEAD_PRONE
	MRI_Default (cryo) / MRI_TxSuc (surface)
Time of start	09:18 AM

TAPPE ears and eyes + do not cover too much the rat and tail not on the heating pad + position well the cryocoll on the region you wish to see -
for MRSI pay attention to STRIATUM

1 Localizer Adj: wobble : Coil 1 VOLUME - setup (manual)-do not do Apply, Coil 2 CRYO: 5Mhz aiwth Start (automatic) - DO PRINT

SCREEN (both)

E2 → Acquire localizer for position with basic frequency/MR scan in the instruction cart
fix the rat

power= **1,61 W**

2 Localizer ref power is adjusted automatically

E3 → Acquire with automatic instructions

Target power for MRSI: 1.6W

3 Loc MS3_Localizer multislice_10_short - for voxel and ellipsoid mapShim

E4 → Seq: 1 avg, 10 slices, Image size: 256x256, FOV: 24x24 → **position the slices**

4 T2 turbo rare 4_T2 turbo rare for voxel position AXIAL (T2 turbo rare_8_52K_250924_Ax)/(T2 turbo rare_6_54K_080822_Sag)

Seq: 2 avg, DS=4, Rare factor: 8, encoding start -1 (in resolution/encoding), 256x256, 28 slices,
0 gap, FOV=24x24

Position the slices, current shim -follow rat temp for all RARE seq

No = **E5** → Automatic acquisition

5 Coronal 5_T2 turbo rare for voxel position CORONAL (T2 turbo rare_8_52K_Cor_mapshim_250924), 20 slices, TE=33ms, TR=2.5s

ACQ the B0 map before in the Adj platform (settings are already saved) - **CEREBELLUM SHOULD BE IN THE FOV**

B0 map No **E7** Can be increased or decreased based on brain size

Mapshim - elliptical 13x9x19 Position the ellipsoid before starting the acquisition

5 Atlas 5_T2 turbo rare_6_52K_Cor_250924_atlas_03_RC - has 2 sat bands left-right -pe direction for aliasing)

Copy slice orientation from one coronal to the other (not the geometry) - **ATTENTION WHEN POSITIONING THE SAT BANDS**

Seq: 10 avg, DS=4, Rare factor: 6, encoding start -1 (in resolution/encoding), 128x128, 40 slices, 0 gap

Position the slices +FOV =24x24 as MRSI

No = **E8** → Automatic acquisition

Create a mapshim report **SD values:**

31 Hz

STEAM

QUALITY CHECK

6_STEAM_water_MRSI_250924 (shim verification following ellipsoid shimming) - "_big" - bigger VOI (+cerebellum)

position voxel - 10 x 2 x 10 mm

Seq: 16 avg, 1 rep, offset 0, no WS, 2 DS, TE=3, TM=10, TR=4

OVS (15/12mm, gap=1mm), spoilers 15-25-35%, ref scan 16, Seq spoilers : 30-36-25%

Current shim

Vox size=10x2x10

position= **0,561x0,477**

x 0,174 mm

No = **E9** → Automatic acquisition

→topspin/tick in Reconstr. **lw = 27,884 Hz**

Save shim (no matter the lw always save the first shim)

6_STEAM_water_MRSI_250924_mapshim - for shimming - copy ONLY voxel geom from previous

If lw is very high (≥25-30Hz) then duplicate 1_Localizer and in Adj platform acquire B0map

B0 map No: **E12**

In STEAM in Setup card: Map shim on the voxel: cuboid, no margin, iterative corrections

No = **E10** → Automatic acquisition (target 15-22 Hz)

Save shim

lw = 27,884 → 17,763 Hz

No = →If necessary duplicate last STEAM acq (if lw is 17-25 Hz): copy seq, Setup/Current Shim, then on Adj platform:
loc.freq/loc.shim/loc.freq - 512 points

Save shim

Final lw = 17,763 Hz

Mandatory: **8_STEAM_met_MRSI_281024_300WS** - copy ONLY voxel geometry from previous

1 repetition with 32 averages, ref scan 16, OVS on (Angles Geometry Display: 0/-180/-180)

WS ON : WS pulses 84/150 → If you change the WS note it here:

Last delay (no 7) in WS - 28ms - test in Setup Mode - **USUALLY IS 28ms, TEST IT BUT ACQ WITH**

THE ONE FROM MRSI - 28.7 (BW=400 Hz) or 29.34 (BW=350Hz)

TE=3ms TM=10ms TR=4s

No = **E13** → Automatic acquisition + Current shim + copy sequence (WS off, 16 av) for H2O acq if absolute quantif is needed

FID-MRSI 9_CSI_fid_cor_metab_TE13_BL4_WScC_100225 - do the traditional ISMRM 2025 sat bands – do a print screen (see end this document) – OLD version: 9_CSI_fid_cor_metab_281024_TE13_BL4_WSOptimized

TR=822ms, novelty of the version WS last delay=28ms (the auto repetition spoiler was decreased to 1.2ms with 35%)

Copy ONLY Slice orientation from Coronal RARE ATALS (scan 5) - check that the values are ok (Geometry tab)

Adjust the slice offset - Middle coordinate of the VOI in STEAM

Adjust the position of saturation bands

Use both Axial and Coronal RARE acq as support

FOV = 24x24 , flip angle = 55, Dummy scans = 8 Linear encoding,

start at -1, -1

weighted acquisition / standard acquisition – Always standard with 1 avg.

Preparation/WS – VAPOR*

Launch in the **setup mode** to find optimal WS and BW of the pulses (test also with 8x8 matrix)

WS flip angles = 84/150

BW = 350-400

Matrix size= 31 x 31

Coronal Slice position (from Atlas):

10 - 16

*For the flip angles 84/150 and BW=400, Hz last delay aprox 28ms works best so start with those

BW=350 Hz can be used if shim is of good quality (17-20Hz for H2O), then last delay will slightly increase to 28.54 ms

Number of averages:

Number of repetitions:

No= E14 → Automatic acquisition + current shim - **metabolites**

Duplicate sequence - preparation/WS - none

No= E15 → Automatic acquisition + current shim - **water**

(Angles Geometry Display: -90/0/-90)

FID-MRSI – 2nd DATA SET (if you do testing for the FOV sat)

!!!! We can copy the previous acquired sequence and just move the sat bands outside the brain (left and right on the coronal, up on the coronal, and obliques one on the coronal)

9_CSI_fid_cor_metab_100225_TE13_BL4_WScC_rooftop - done on 100225 – lipid suppression is better and we can keep the sat bands outside the brain

WS -last delay=29.4ms for 350HZ, WS -last delay=29.3ms for 400HZ

If you repeat with other parameters

High Res

1st repeat

FID - MRSI

47 x 47

weighted acquisition / standard acquisition - circle the one you select

Preparation/WS – VAPOR

Averages: 3

Repetitions: 1

No= E16 → Automatic acquisition + current shim - **metabolites**

No= E17 → Automatic acquisition + current shim - **water**

CS: AF=4 / Core=20%

2nd repeat

FID - MRSI

63 x 63

weighted acquisition / standard acquisition - circle the one you select

Preparation/WS - VAPOR

Averages = 3

Repetitions= 1

No= E18 → Automatic acquisition + current shim - **metabolites**

No= E19 → Automatic acquisition + current shim - **water**

CS: AF=4 / Core=20%

CENTER THE MRSI SLICE ON THE CORONAL IMAGE

E20 - STEAM for 3D Shim : $hw = 24 \text{ Hz}$

3D MRSI -- checked on 04/02/2025 -OK

13_cristinaCSI 3D 100225 FOV BL1 WSc - WS=400Hz, last delay 29.13

SHIM THE 3D VOLUME first

Copy slice orientation from Atlas (not 2D MRSI), check in Routine that slice $thk=1\text{mm}$, slice offset – take the one from 2D+2mm to have the 3rd slice in the middle of the 2D slice

$Met = 21$

$Wst = 22$

2D Multislice

Multislice

13_cristinaCSI Multislice 100225 FOV BL1 as3DWS

SHIM THE 3D VOLUME first

- Copy sat bands from 3D, Routine slice $thk=1\text{mm}$, copy slice orientation and geometry from 3D, again check thk slice in Routine, In geometry select 2D, In geometry put same slice offset as in 3D and keep 1 slice package but scroll to add 9 slices
- $7398 \text{ ms} = TR$, 90RF, 8 DS, $WS=400\text{Hz}$, last delay 29.5. ms, Pay attention to have the same WS param as in the previous 2D scans – same BW
- You need to import sat bands from 3D
- Auto rep spoiler 1ms, 50%

$Met = 23$

$Wst = 24$

ATTENTION – DURATION OF AUTO REPETITION SPOILER (SEQUENCE TAB) INCREASES THE LAST DELAY IN WS

MM acquisition protocol -- checked on 04/02/2025 -OK

13_cristinaCSI_DIR_MM_100225_RFandWSc_BA_BL_cc -done in 100225 to keep the smallest last delay in WS (28ms for 400Hz BW)

BA_djCSI_excpul_DIR_MM_100225_90RF_TE06_BL4_WSc --- $AD=0.6\text{ms}$, $TR=3600\text{ms}$, Gav-standard, 2048p, 5000BW, $WS=400\text{Hz}$ – 28ms, P10- 0.6ms for 90RF,

ATTENTION !!!

- WS bandwidth will range between 350 – 440 Hz depending on the quality of the shimming
- Duration of the TI for the DIR modules is $TI1=2200 \text{ ms}$ and $TI2= 850$ or 750 ms

Short AD protocol --- checked on 04/02/2025 -OK

10_djCSI_excpul_100225_06_BL4_WSc -P10 (0.4ms due to power req), $pe=0.2$

10_CSIfid_cor_metab_100225_TE065_BL4_WSc – RF=0.3ms, $pe=0.2\text{ms}$

PRESS MRSI – checked on 04/02/2025 -OK

9_CSI_PRESS_cor_metab_100225

$TE=10.2\text{ms}$ as 14T, $TR=2\text{sec}$; 32 min

Same Spectroscopy tab as FID MRSI (768p, 5000Hz BW)

$10 \times 10 \times 2$ – (-90, 0, -90) - for FID-MRSI (-90,0-90)

WS as FID-MRSI just that last delay is 28ms for 350Hz – should we use the same as for FID-MRSI?? (if yes then we lose the only advantage of this seq)

OVS as for STEAM, 1mm gap, 10mm slab, 3ms spoilers (15-25-35%), RF 1ms sech

Seq tab: 90 RF – 0.5ms (8400Hz) but 180 RF -0.6ms due to power demands (4250Hz)

-2ppm as FID-MRSI; auto Rep Spoilers (2ms, 25%), Auto Echo SPolire (1.5ms, 25%) – with this min $TE=9.8\text{ms}$ –

is this spoiler is decreased then TE is also decreased ☺

For 3D & Multislice :

slice 1 = 6

slice 2 = 10

slice 3 = 13

slice 4 = 16

slice 5 = 20

slice 6 = 23

slice 7 = 26

slice 8 = 30

slice 9 = 33

Application ---- MRS SVS cryo

STEAM QUALITY CHECK and SHIM

1. **8_2_STEAM_water_hippo_281024_OVS** - position the voxel in the hippo and check the water linewidth after the shims done before (**usually we get between 14-18Hz**)
2.8x2x2.8 - position=
H2O - lw =
WHEN SCANNING SICK ANIMALS THIS SCAN IS NOT NEEDED
2. **8_2_STEAM_water_hippo_281024_OVS_mapshim** - copy ONLY Voxel geometry - shim the voxel to improve (**target values: 9-12 Hz**)
H2O lw =
Save shim
3. **Duplicate the previous sequence if lw is ≥ 11 Hz** - in the Setup Card select Current shim - in the adj platform do 1) Local freq; 2) local Shim; 3) Local freq by selecting in the Spectroscopy card only 512 points in the FID
Acq the signal
H2O lw =
Save shim

WS CALIBRATION

4. **9_STEAM_met_hippo_281024_OVS_32av** - test WS (**BW=270 Hz and last delay=28ms**) - copy ONLY voxel geometry from previous scan
In Setup Mode you can test the WS
WHEN SCANNING SICK ANIMALS - JUST USE SETUP MODE TO TEST THE WS YOU DO NOT NEED TO ACQ THE DATA, however you need to copy the voxel and WS param if changed in the next scan

ACQ OF METABOLITE AND 2XWATER SIGNAL

5. **9_STEAM_met_hippo_281024_OVS_16x16** - acq metabolite signal - copy ONLY voxel geometry from previous scan
Check that the voxel position is the good one
6. **9_STEAM_met_water_hippo_281024_OVS_absQ** - acq water signal for absolute quantif with OVS
Copy ONLY voxel geometry from previous scan
Check that the voxel position is the good one
7. **Duplicate -- 9_STEAM_met_water_hippo_281024_OVS_absQ** - take out OVS and acq water signal for absolute quantif without OVS

sSPECIAL

Specialjm_adiabIR_hippo_metf_281024

Specialjm_adiabIR_hippo_water_absQ_281024

Sat Bands - the strict version - \\cibmaitsrv1.epfl.ch\scannerdata\Toi\20250113